



The release notes describe what’s new for the Logistics Portal. In version 2022.01 of the Logistics Portal, the following updates are available:

Medioh API Integration.....	1
Glucometer Term Remediation.....	2
Update Kit Sync Report – Highlight unmatched columns	5
Kit Type Master List Updates.....	5
Record Information Card – Cuff Size and Scale Type	7
Fixes.....	7

Medioh API Integration

Several APIs were updated or created to support integrating the Logistics Portal and Medioh via APIs for fulfilling kit shipping, kit pickup, kit reprocessing, device orders, and device replacement orders.

New APIs:

- **Fulfillment API for Kit Ship** (Outbound) – Fulfillment API was created for Medioh so that Kit Ship fulfillment related data can be sent to Medioh when an order is created in the Logistics Portal.
- **Pick Slip Drop API** (Inbound) – Pick Slip Drop API was created so that kit inventory and kit ship status can be received by Logistics Portal from Medioh during preparing the kit for shipment. Upon receipt of this API information, Kit Ship order status is marked as In Progress.
- **Device Linking API** (Inbound) – Device linking API was created so that Medioh can send the kit number and assigned devices/parts and their Mac ID (Bluetooth address) to be linked to the fulfillment order.

- **Ship Confirmation API for Kit Ship orders** (Inbound) – Kit Ship Confirmation API was created so that kit ship information can be received by the Logistics Portal from Medioh when a kit is shipped.
- **Pickup Order Alert API** (Outbound) – Pickup Order Alert API was created so that Medioh is alerted to receive the incoming kit/device/part. This API is triggered when a pickup order is created.
- **Cancel Order API** (Outbound) – Cancel Request API was created so that related data can be sent to Medioh when a kit ship order is canceled in the Logistics Portal.

Updated APIs:

- **Fulfillment API for Kit Ship** (Outbound) – The following additional fields were added to the Fulfillment API: Selected BP Cuff Size, Selected Scale Type, Caregiver Name, Caregiver Phone, and Caregiver Email
- **Kit Ship Confirmation API** (Inbound) – Ship Confirmation API has the following updates:
 - Products included in the call are listed in the array.
 - Device SKU/part number for the SIM can be included in the product field.
- **Device Linking API** (Inbound) – Device Linking API has the following updates:
 - Support for non-kits on demand orders when a kit is swapped.
 - Device SKU/part number for the SIM can be included in the product field.
 - Any changed devices or parts are removed individually instead of removing them all.
 - Products included in the call are listed in the array.



Glucometer Term Remediation

Due to the word *Glucometer* being a trademarked term, several instances of this word in the Logistics Portal were changed to *Glucose Meter*.

Kit List Page


On the Kit List page, the following columns were updated:

- Glucometer to Glucose_Meter
- Glucometer_Part_Num to Glucose_Meter_Part_Num
- Glucometer_Serial to Glucose_Meter_Serial
- Glucometer_Bt_Address to Glucose_Meter_Bt_Address




Kit Details Page

On the Kit Details page, **Glucometer** was renamed to **Glucose Meter**.




Devices Page

On the Devices page, **Glucometer** was updated to **GlucoseMeter**.




Reprocessing Report

On the Reprocessing report, the **Glucometer** column was updated to **Glucose Meter**.



AT&T Kit Inventory Report


On the AT&T Kit Inventory report, the column **Glucometer** was renamed to **Glucose Meter**.



Update Kit Sync Report – Highlight unmatched columns


In the Kit Sync report, columns that are out of sync/unmatched for kits will be highlighted to see which categories/properties are unmatched quickly.

Note: If there are no differences, then no highlighting is applied.




Kit Type Master List Updates

On the Kit Type Master List, the Kit Type Names and Descriptions were updated. Some kit type names that contained **Glucometer** were updated to **Glucose Meter**.



On the **Kit Type Master List > Add/Edit Kit Type** window, a new **Icon** field was added. The **Icon** field is optional and read-only.

Note: The process of syncing the Kit Type data from the Logistics Portal to the Care Team Portal was updated to include the new Kit Type Icons.



Record Information Card – Cuff Size and Scale Type

On the Fulfillment Order Details page, **BP Cuff Size** and **Scale Type** were added.

- If the Kit Type for the fulfillment order includes a blood pressure monitor device, the selected BP cuff size will appear in the **Fulfillment Order > Record Information**.
- If the Kit Type for the fulfillment order includes a scale, the selected scale type will appear in the **Fulfillment Order > Record information**.

Record Information

Fulfillment Type: Kit Ship
Kit Pool: Default Kit Pool
Kit Type: Tablet + BP + PulseOx + Glucose Meter
Kit Id: 100000
Logistics Center: Vivify Development
BP Cuff Size: Wide Range
Scale Type: Normal

Contact Attempts: 0

Update Contact Attempts

Agent:
Assign Agent

Fixes

- **VIVY-6758** – Logging was added to KitShippingController and KitShippingComponent when creating a Kit Ship Request in order to provide meaningful data when an error requires diagnosis.
- **VIVY-7648** – When the Update Contact Attempts button was clicked on a Kit Pickup detail page, a Fulfillment Note was saved, but the Contact Attempt Date and Note were not saved to the kit shipping info database. This functionality was updated to save as expected.
- **VIVY-7682** – The Kit detail page showed a blank location for the current kit location in some scenarios. This occurred when the kit did not have a history of being received. The page was updated to display the current kit location.
- **VIVY-7519** – Performance improvements were made for Add Device API.
- **VIVY-6678** – The blue cards on Logistics > Dashboard page for the Kit Ship and Kit Device Replacement columns were not visible. The Dashboard was updated to display these cards.
- **VIVY-6679** – An error was appearing when searching on the Kit Storage page. This error was resolved, and the search is now functioning as expected.

- **VIVY-7507** – The Ship API was updated to return an error when a re-processing record exists.
- **VIVY-7629** – Errors that appeared in the create database process for test portals were resolved.
- **VIVY-7361** – When a user had all logistics locations access, this caused an empty list to appear. The system was updated to show the complete list on the PO page in this situation.
- **VIVY-7728** – When deleting a user from the User View page, a confirmation or warning did not appear. The confirmation window was added.
- **VIVY-5199** – A migration issue that occurred during deployment was resolved.
- **VIVY-5046** – When creating a user in the Logistics Portal, a space before the username, first name, or last name caused an error. Validation was added to prevent the error.
- **VIVY-3147** – When shipping a kit, the Name field was not required. The Name field is now required for shipping.