



Release Notes

Version 2021.03 of the Logistics Portal includes system updates to support the Logistics team and their work in the portal.

System Updates

System Update	Description
VIVY-5353	The Mass Kit Sync/Update feature was not sending the correct status ID when syncing kits. For example, when syncing a kit with the Stored status, the system was sending an On Hold status. The system was updated to sync statuses as expected.
VIVY-5313	The Contact Attempt message and count was not syncing to the Care Team Portal from the Logistics Portal. The portal was updated to sync Contact Attempts.
VIVY-2914	The Create Kits page was updated to use Vue, so it is in line with the other new pages in the portal and the functionality is improved. The page will function the same with minor updates. For example, adding parts is now a single click and the page will now enforce that the user only creates the number of kits that are on the order.
VIVY-4920	On the Kit List page, when updating a group of kits to a new status for a second time, it wasn't possible to complete the action without refreshing the page. The page was updated to load as expected when a second update is completed.