



Release Notes Logistics™


December 2021

The release notes describe what's new for the Logistics Portal. In version 2021.12 of the Logistics Portal, part numbers are visible and additional logging was added for APIs.

Part Numbers

Part numbers are now visible in the Logistics Portal as required by the Food and Drug Administration (FDA). The part number displays to the left of the version number in the footer. The part number can be used as an additional identifier when contacting Vivify Support, if needed.

Example: Logistics Dashboard with Part Number



Fixes

- **VIVY-6758** – Logging was added to KitShippingController and KitShippingComponent when creating a Kit Ship Request. This update will allow the API to return more meaningful messages, which is helpful for troubleshooting when an error occurs.
- **VIVY-7607** – Additional logging was added for the Device Receipt API to assist with troubleshooting issues.
- **VIVY-7632** – An error message appeared when swapping kit devices with the Ship Confirm API. The API was updated to resolve this error.